MiniLink[®] REAL-Time

TRANSMITTER

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TRANSMITTER

TRANSMISOR

ZENDER

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TRANSMITTER

LÄHETIN

SÄNDARE

SENDER

TRANSMISSOR

ODDAJNIK

VYSÍLAČ

TÁVADÓ

NADAJNIK

VERİCİ

VYSIELAČ

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Azerbaijan: Isomed
Tel: +994 (12) 464 11 30

Bangladesh: Sonargaon Healthcare Pvt Ltd.
Mobile: (+91)-9903995417
or (+880)-1714217131

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Japan: Medtronic Japan Co. Ltd.

Tel: +81-3-6430-2019
24 Hr. Support Line: 0120-56-32-56

Kazakhstan: Medtronic BV in Kazakhstan

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Malaysia: Medtronic International Ltd.

Tel: +603 7946 9000

Middle East and North Africa: Regional Office

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Sri Lanka: Swiss Biogenics Ltd.
Mobile: (+91)-9003077499
or (+94)-777256760

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
RELEASED

RELEASED

The Medtronic MiniLink® transmitter is a component of select continuous glucose sensing systems and sensor-enabled insulin pump systems. The transmitter powers the glucose sensor, collects glucose data, and wirelessly sends the data to a Medtronic continuous glucose sensing monitor or Medtronic MiniMed® sensor-enabled insulin pump.

A complete MiniLink kit includes:

- MiniLink transmitter (MMT-7703)
- Watertight tester (MMT-7726)
- Sensor insertion device
- Charger (MMT-7715)
- AAA or LR-03 alkaline battery(ies)

Indications for use

The transmitter is indicated for use as a component of select Medtronic continuous glucose sensing systems and Medtronic MiniMed sensor-enabled pump systems.

Contraindications

None known.

Warnings

Product contains small parts and may pose a choking hazard for young children.

The sensor should be removed if redness, bleeding, pain, tenderness, irritation, or inflammation develops at the insertion site, or if you experience unexplained fever.

Optional occlusive dressing should be removed if irritation or reaction to the tape develops.

Wait five minutes after sensor insertion before connecting the transmitter to the sensor.

- Make sure that the site is not bleeding before connection.
- If bleeding occurs, apply steady pressure with a sterile gauze or clean cloth at the insertion site until bleeding stops. After bleeding stops, attach the transmitter to the sensor.

- If bleeding persists after three minutes, remove the sensor and discard. Insert a new sensor in a different location.

Contact your local help line or representative for assistance if you experience any adverse reactions associated with the transmitter or sensor.

Exposure to magnetic fields and radiation

If you are going to have an X-ray, MRI, diathermy treatment, CT scan, or other type of exposure to radiation, take off your sensor and transmitter before entering a room containing any of these equipment. The magnetic fields and radiation in the immediate vicinity of these devices can make the transmitter non-functional. If your transmitter is inadvertently exposed to a magnetic field, discontinue use and contact your local help line or representative for further assistance.

Airport security

Important information about airport security systems, and using your transmitter on an airplane, can be found on the Emergency Card. Be sure to carry the Emergency Card provided when you are traveling.

Precautions

Establish a rotation schedule for choosing new sensor sites. Avoid sites that are constrained by clothing, have scar tissue, or are subject to rigorous movement during exercise.

Notice

This device complies with the United States Federal Communications Commission (FCC) and international standards for electromagnetic compatibility.

This device complies with Part 15 Rules. Operation is subject to the following two conditions: (1) this device may not cause harmful interference, and (2) this device must accept any interference received, including interference that may cause undesirable operation.

The transmitter does not interfere with any radio frequency signals transmitted from outside sources. These FCC standards are designed to provide reasonable protection against excessive radio frequency interference and prevent undesirable operation of the device from unwanted electromagnetic interference.

Important: Changes or modifications not expressly approved by the party responsible for compliance could void the user's authority to operate the equipment.

RF interference from other devices

Common consumer electronic devices that transmit in the same frequency band used by the transmitter may prevent the receiving device (Guardian® REAL-Time monitor or MiniMed Paradigm® REAL-Time insulin pump) from receiving the glucose information sent by the transmitter. Most cellular (mobile) phones and 900 MHz cordless phones, when transmitting or receiving, may cause significant interruption of transmitter-receiver communication. It is likely that other devices operating in similar frequency ranges will have a similar effect. This interference, however, will not cause any incorrect data to be sent and will not cause any harm to your transmitter.

Guardian REAL-Time monitors and MiniMed Paradigm REAL-Time insulin pumps include a programmable “Weak Signal” alert that notifies you when one or more transmitter transmissions were not received as expected by the receiving device. (The receiving device will also issue a “Lost Sensor” alert if communication is interrupted for approximately 40 minutes.)

Communication problems can typically be resolved by ensuring that the distance between transmitter and receiving device is less than 1.8 meters (six feet), and by turning off or moving away from other RF transmitting devices. You can also reorient or relocate the transmitter and/or the receiving device to try to correct the interference. Testing conducted with several different cellular phones suggests that interference will not be a problem if the phone is at least 31 cm (12 inches) from the transmitter or receiving device while it is being used (greater separation distance may be required for certain devices).

Medtronic declares that this product is in conformity with the essential requirements of Directive 1999/5/EC on Radio and Telecommunications Terminal Equipment.

For additional information, contact Medtronic MiniMed at the address or phone number provided on the back cover.

Assistance

Please contact your local representative for assistance. Refer to the Medtronic Diabetes International Contacts list at the beginning of this user guide for contact information.


Charger

The transmitter contains a non-replaceable, rechargeable battery that you can recharge as needed with the charger. The charger has a green light that shows the charging status and a red light that communicates any problems during charging. If you see a red light, see the Troubleshooting section. The charger needs one AAA alkaline battery, size E92, type LR03 (Energizer brand recommended) to operate.


Note: A new AAA or LR-03 battery contains enough power to recharge the transmitter more than 40 times. If the battery is installed incorrectly or is low, the charger will not work. Repeat the battery installation steps using a new battery.

Installing a battery in the charger

- 1 Push the battery cover in and slide it off (as shown).

- 2 Insert a new AAA or LR-03 battery. Make sure that the + and - symbols on the battery align with these same symbols shown on the charger.
- 3 Slide the cover back on the charger until it clicks into place (as shown).



Charging the transmitter

Before using the transmitter for the first time, you must fully charge the transmitter battery. This may take up to eight hours. It is recommended to recharge the transmitter after each sensor use. If you choose to recharge the MiniLink after a 6-day sensor use,

the charging time will be less than 20 minutes. A fully charged transmitter battery works more than 14 days without recharging. After 14 days of use, the transmitter recharges in less than two hours.

To charge the transmitter:

- 1 If a green light on the transmitter is lit or flashing, do not connect it to the charger. The transmitter will not charge with its green light on. Wait until the green light turns off (approximately 30 seconds), and then connect the transmitter to the charger.
- 2 Connect the transmitter to the charger by lining it up, flat side down, with the charger. Push the two components together fully.
- 3 Within 10 seconds after the transmitter is connected, a green light on the charger will flash for 1 to 2 seconds as the charger powers on. For the rest of the charging time, the charger's green light will continue to flash in a pattern of four flashes with a pause between the four flashes.
- 4 When charging is complete, the green light on the charger will stay on, without flashing, for 15 to 20 seconds and then turn off.
- 5 After the green charger light turns off, disconnect the transmitter from the charger. The green light on the transmitter will flash for about 5 seconds and then turn off.

Setting up the transmitter

For instructions on setting up your transmitter, see the user guide for your pump.

Connecting the transmitter to the sensor

After you have confirmed that the insertion site is not bleeding, connect the transmitter to the sensor.


Note: *The transmitter must be fully charged before connecting it to the sensor.*

Before physically connecting your transmitter to your sensor, consult the user guide for your pump or monitor for details on setting up and starting the sensor.

Once connected, the sensor requires a short warm-up period. The pump notifies you when the warm-up period ends and the sensor is ready to use.


To connect the transmitter to the sensor:

- 1 If you are using Enlite Sensor (MMT-7008), apply the Enlite overtape to the sensor and wait five minutes before connecting the transmitter. Consult the user guide for your Enlite Serter for details on applying Enlite overtape. If you are using a different sensor (MMT-7002/MMT-7003), wait 10 to 15 minutes before connecting the transmitter.
- 2 Touch the rounded end of the inserted sensor to prevent it from moving during connection.
- 3 Hold the transmitter as shown. Line up the two notches on the transmitter with the side arms of the sensor. The flat side of the transmitter should face the skin.
- 4 Slide the transmitter onto the sensor until the sensor's flexible arms snap into the notches on the transmitter. If the transmitter is properly connected, and if the sensor has had enough time to become hydrated, the green light on the transmitter will flash within 10 seconds.
- 5 If the transmitter light does not flash, disconnect it from the sensor, wait for several seconds and then reconnect. If the transmitter light still does not flash, charge the transmitter.
- 6 After the transmitter light flashes green, use your pump or monitor to communicate with the sensor. For more instructions, see the user guide for your pump or monitor.
- 7 After the transmitter successfully sends sensor data to the pump or monitor, attach the adhesive tab from the Enlite sensor (MMT-7008) to the transmitter. For other sensors (MMT-7002 or MMT-7003), continue to the next step.
- 8 [Optional]: Apply occlusive dressing over the transmitter and the sensor.

Disconnecting the transmitter from the sensor

- 1 Carefully remove any occlusive dressing from the transmitter and sensor.
- 2 For the Enlite sensor, remove the adhesive tab from the top of the transmitter.
- 3 Hold the transmitter as shown, and pinch the flexible side arms of the sensor between your thumb and forefinger.
- 4 Gently pull the transmitter away from the sensor.

Bathing and swimming


After the transmitter and sensor are connected, they form a waterproof seal to a depth of 2.4 meters (eight feet) for up to 30 minutes. You can shower and swim without removing them. No occlusive dressing is needed.

Watertight tester

The tester is used to test the transmitter to make sure it is working. It is also used as a required component for cleaning the transmitter. Properly connecting the tester to the transmitter will ensure that fluids do not come in contact with the transmitter's connector pins. Fluids can cause connector pins to corrode and affect the transmitter's performance.


Each tester can be used a total of 30 times. Open and use one tester at a time. Keep track of the tester uses on the supplied log sheet and discard the tester after 30 uses in a bio-waste container. If you continue to use the tester beyond 30 times, the transmitter's connector pins could be damaged, because the tester cannot continue to provide a waterproof seal. For instructions on how to check the connector pins, see *Inspecting the transmitter connector pins*, on page 7.

Caution: Only use the watertight tester with the transmitter. Do not use any other test plug.

Inspecting the transmitter connector pins

This image is an example of how the connector pins should look.


Look inside the transmitter's connector opening to make sure that the connector pins are not damaged or corroded. If the connector pins are damaged or corroded, the transmitter cannot communicate with the charger or pump. Contact your local help line or representative for assistance. It may be time to replace your transmitter.

Also look for moisture inside the connector opening. If you see any moisture, allow the transmitter to dry for at least one hour. Moisture inside the connector opening could cause the transmitter to not work properly, and could cause corrosion and damage over time.

Connecting the tester for testing or cleaning

To connect the tester:

- 1 Hold the transmitter and the tester as shown. Line up the flat side of the tester with the flat side of the transmitter.
- 2 Push the tester into the transmitter until the flexible side arms of the tester click into the notches on both sides of the transmitter. Within 20 seconds the green light on the transmitter will flash for about 10 seconds when properly connected.
- 3 To test the transmitter, check the sensor icon on the pump or monitor to ensure that the transmitter is sending a signal (see the user guide for your device).
- 4 To clean the transmitter, see *Cleaning/disinfecting the transmitter*, on page 9.
- 5 After testing or cleaning, disconnect the tester from the transmitter.


Disconnecting the tester

To disconnect the tester:

- 1 Hold the transmitter body as shown and pinch the side arms of the tester.
- 2 With the tester arms pinched, gently pull the transmitter away from the tester.

Note: To save transmitter battery life, do NOT leave the tester connected after cleaning or testing.

Cleaning/disinfecting the transmitter

The transmitter is intended for personal use at home (single-patient use) or for use in healthcare facilities (multiple-patient use). Single-patient use requires cleaning after each use, while multi-patient use requires cleaning and disinfection after each use. When using the transmitter in a healthcare facility, always follow the cleaning and disinfecting procedure for multiple-patient use.

Caution: Do not discard the transmitter in a medical waste container or otherwise subject it to incineration. The transmitter contains a battery that may explode upon incineration.

Note: The tester is a required component for cleaning and disinfecting the transmitter. For details, see *Watertight tester*, on page 7.

For single-patient use

Always clean the transmitter after each use.

To clean the transmitter:

- 1 Wash your hands thoroughly.
- 2 Attach the tester to the transmitter.
- 3 Rinse the transmitter under cool tap water for at least one minute and until visibly clean. Make sure all hard-to-reach areas are rinsed completely.
- 4 Prepare a mild liquid soap solution using 5 milliliters (one teaspoon) of mild liquid soap per 3.8 liters (one gallon) of room temperature tap water. Never use organic solvents, such as paint thinner or acetone, to clean the transmitter.
- 5 With tester still attached, submerge the transmitter in the mild liquid soap solution and soak for one minute.
- 6 Holding the tester, brush the entire surface of the transmitter using a soft-bristled toddler toothbrush. Make sure to brush all hard-to-reach areas until visibly clean.
- 7 Rinse the transmitter under running room temperature tap water for at least one minute, and until all visible liquid soap is gone.
- 8 Dry the transmitter and tester with a clean, dry cloth.
- 9 Place the transmitter and tester on a clean, dry cloth and air dry them completely.
- 10 Disconnect the tester from the transmitter.

For multiple-patient use

When using the transmitter in a healthcare facility, always clean and disinfect the transmitter after each use.

To clean and disinfect the transmitter:

- 1 Wash your hands and put on gloves.
- 2 Inspect the inside of the transmitter's connector opening for any sign of body fluid. For instructions on how to inspect the connector pins, see *Inspecting the transmitter connector pins*, on page 7.

Caution: The person inspecting the transmitter must have sufficient vision that enables him or her to see small drops of body fluid or debris.

Warning: If you see any body fluid in the connector opening, you must discard the transmitter. Because the transmitter contains a battery, do not discard in a bio-waste container. Instead, continue to clean and disinfect the transmitter, and then discard according to local regulations for battery disposal (non-incineration).

- 3 Attach the tester to the transmitter.
- 4 If there is adhesive residue on the transmitter, follow the instructions in *Removing adhesive residue*, on page 11.
- 5 Rinse the transmitter under cool tap water for at least one minute and until any visible debris is gone.
- 6 Prepare an enzymatic solution using 29.5 milliliters (1 ounce) of enzymatic detergent per 3.8 liters (1 gallon) of lukewarm tap water. Make sure to prepare a fresh solution for each use.
- 7 With the tester still attached, fully submerge the transmitter in the enzymatic solution for one minute.
- 8 Holding the tester, remove the transmitter from the solution. Brush the entire surface of the transmitter using a soft-bristled brush, paying close attention to hard-to-clean areas, until visibly clean.
- 9 Rinse the transmitter under room temperature tap water until all visible detergent is gone.
- 10 Dry any excess moisture by wiping the outside of the transmitter with a clean, dry cloth.
- 11 Prepare a 1:10 bleach solution by using one (1) part 8.25% bleach to nine (9) parts water, for a final concentration of 0.825%. Make sure to prepare a fresh solution for each use.

- 12 With the tester still attached, soak the transmitter in the bleach solution for 20 minutes.
- 13 Rinse the transmitter under room temperature tap water for three minutes.
- 14 Holding the tester, wipe the transmitter with 70% isopropyl alcohol.
- 15 Place the transmitter and tester on a clean, dry cloth and air dry them completely.
- 16 Important: If you saw any body fluid inside the connector opening on earlier inspection, you must now discard the transmitter with tester still attached, according to local regulations for battery disposal (non-incineration).
- 17 Disconnect the tester from the transmitter by gently squeezing the arms of the tester.
- 18 Inspect the housing of the transmitter for any signs of cracking, flaking, or damage. If you see any of these signs, you must now discard the disinfected transmitter according to local regulations for battery disposal (non-incineration).

Warning: Cracking, flaking, or damage of the housing are signs of deterioration and the performance of the device may be compromised. This may affect the ability to properly clean and disinfect the transmitter. If these signs are noted, the device must be discarded according to local regulations for battery disposal (non-incineration).

- 19 Discard the used gloves and thoroughly wash hands with soap and water.

Removing adhesive residue

You may need to perform this procedure only if you have used optional occlusive dressing, which may leave adhesive residue on the transmitter. If you visually inspect the transmitter and see adhesive residue on it, follow the instructions below.

To remove adhesive residue, you will need the following materials: Detachol® medical adhesive remover and cotton swabs. You can buy Detachol at <http://www.amazon.com/>.

To remove adhesive residue:

- 1 Make sure the tester is attached to the transmitter.
- 2 Holding the tester, saturate a cotton swab in the Detachol solution and gently rub the adhesive residue on the transmitter until it is fully removed.
- 3 Continue with the appropriate single-patient or multi-patient cleaning procedure above.

Cleaning the charger

The charger cannot be disinfected. This procedure is for general cleaning as required, based on physical appearance.

Caution: The charger is **NOT** waterproof. Do **NOT** immerse in water or any other cleaning agent.

Caution: Dispose the charger according to the local regulations for battery disposal (non-incineration).


Warning: For multiple-patient use, always clean and disinfect the transmitter after removing it from the patient and before attaching it to the charger. If blood comes in contact with any surface areas of the charger, the contaminated device must be discarded. The charger contains a battery which may explode upon incineration.

To clean the charger:

- 1 Wash your hands thoroughly.
- 2 Use a damp cloth with mild cleaning solution, such as a dishwashing detergent, to clean any dirt or foreign material from the outside of the charger. Never use organic solvents, such as paint thinner or acetone, to clean the charger.
- 3 Place the charger on a clean, dry cloth and air dry for 2–3 minutes.


Troubleshooting

Question: Why do I see quick flashing red lights on the charger?


Answer: The transmitter battery is very low. Leave the transmitter on the charger for 8 hours to completely recharge. If the red light is still flashing after 8 hours, leave your transmitter on the charger for 24 hours. If the red light is still flashing after it has been charging continuously for 24 hours, contact your local help line or representative for assistance. It may be time to replace your transmitter.

Question: Why did the flashing green charger light turn off and a longer flashing red charger light turn on during charging?

Answer: Your charger battery is low. Make sure the transmitter is not connected to the charger and replace the charger battery with a new AAA or LR-03 battery.

Question: Why do I see a mix of quick and long flashing red lights on the charger?

Answer: Your charger AND transmitter batteries are very low. Replace the charger's AAA or LR-03 battery. If you now get the pattern for very low transmitter battery, leave the transmitter on the charger for 8 hours to recharge. If the red light is still flashing after 8 hours, leave your transmitter on the charger for 24 hours. If the red light is still flashing after it has been charging continuously for 24 hours, contact your local help line or representative for assistance. It may be time to replace your transmitter.

Question: I had my transmitter on the charger for a day. Will this damage my transmitter?

Answer: It will not damage the transmitter. You cannot overcharge it.

Question: What should I do if the transmitter green light does not flash when connected to the sensor?

Answer: Is the sensor inserted in the body? If **it is not inserted**, the transmitter will not flash green or send signals to your pump or monitor.

If the **sensor is inserted** in the body, you need to disconnect the transmitter from the sensor, wait for several seconds and then reconnect. If the green light still does not flash, charge the transmitter.

Question: Why didn't I see the transmitter green light flash after connecting it to the tester?

Answer: Check the connection. If you still do not see a green light flash, fully recharge the transmitter battery. Test the transmitter with the tester. If you still do not see a green light flash, contact your local help line or representative for assistance. It may be time to replace your transmitter.

Storing the devices

Store the transmitter, charger, and tester in a clean, dry location at room temperature. Although not required, you may store the transmitter on the charger. If the transmitter is not in use, you must charge the transmitter at least once every 60 days.

Specifications

Biocompatibility	Transmitter: Complies with ISO 10993-1 for body contact
Applied parts	Transmitter Sensor
Operating conditions	Transmitter Temperature: +32° to +122°F (0° to 50°C) Caution: When operating the transmitter on a tester in air temperatures greater than 41°C (106°F), the temperature of the transmitter may exceed 43°C (109°F). Transmitter Relative Humidity: 10% to 95% with no condensation Charger Temperature: +50° to +104°F (10° to 40°C) Charger Relative Humidity: 30% to 75% with no condensation
Storage conditions	Transmitter Temperature: -4° to +131°F (-20° to +55°C) Transmitter Relative Humidity: 10% to 100% with condensation Charger Temperature: +14° to +122°F (-10° to +50°C) Charger Relative Humidity: 10% to 95% with no condensation
Battery life	Transmitter: 14 days of continuous glucose monitoring immediately following a full charge Charger: Completes 40 typical charging operations using a new AAA or LR-03 battery
Transmitter frequency	MMT-7703NA 916.5 megahertz MMT-7703WW 868.35 megahertz
Declaration of Conformity	Medtronic declares that this product is in conformity with the essential requirements of Directive 1999/5/EC on Radio and Telecommunications Terminal Equipment and Directive 93/42/EEC on Medical Devices.

Icon Table

	Serial number
	Model number
(1X)	One per container/package
(2X)	Two per container/package
	Date of manufacture
	Manufacturer
	Refer to instruction manual/booklet before every use. (Appears blue on label.)
	[Storage] Temperature limit
	Marking of [European] conformity. This symbol means that the device fully complies with MDD 93/42/EEC (NB 0459) and R&TTE Directive 1999/5/EC.
	RF device distributed in Australia
	Non-ionizing electromagnetic radiation
CONF	Configuration
	Type BF applied part: Degree of protection against electric shock
IC	RF device distributed in Canada
IPX8	Degree of protection against water under conditions identified by the manufacturer: 8 feet or 2.4 meters immersion for 30 minutes or less.
	Caution: Operator awareness is required in order to avoid undesirable consequences.
	[Storage] Humidity limitation
	Authorized representative in the European community

©2014 Medtronic MiniMed, Inc. All rights reserved.

MiniMed® and MiniLink® are registered trademarks of Medtronic MiniMed, Inc. Enlite™ is a trademark of Medtronic MiniMed, Inc. Detachol® is a registered trademark of Ferndale Laboratories Inc. Energizer® is a registered trademark of Eveready Battery Company.

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Medtronic B.V.